

SION MEDICAL Ltd.

Advancing Healthcare Through Innovation

Products Catalogue

Skin Cleansing
& Disinfection



Lubricating
Jelly



Skin Treatment



Advanced
Wound Care



Bandages, Pads
& Absorbent
Dressing




Packs, Sets
& Kits



General Supply



Medical Devices



PPE
Personal Protective Equipment



January 2024

Sion Biotext Medical Ltd.
www.sn-medical.com

3 Eshel St. Sapirim Industrial Park
P.O. Box 109, 8701002 Sderot, Israel

Tel: +972-8-6898650
Fax: +972-8-6898572

sales@sn-medical.com
customer-service@sn-medical.com

Table of Contents

Skin Cleansing & Disinfection3-7

Lubricating Jelly8-9

Skin Treatment10-12

Advanced Wound Care13-16

Bandages, Pads & Absorbent Dressing17-24

Packs, Sets & Kits25-31

General Supply32-33

Medical Devices - Nilymed Line34-40

PPE - Personal Protective Equipment41-42

About Us

Sion Medical Ltd. is a leading manufacturer of single-use skin cleansing and dermatologic treatment, advanced wound care, burns remedy and surgical solutions for the professional healthcare sector.

Sion Medical Ltd. currently operates through 6 locations worldwide, 5 state-of-the-art production facilities, as well as two sales and logistic companies in Europe and in North America.

The company has 4 decades of experience in the production of Medical Devices and OTC Drugs. Sion Medical Ltd. is on the front line of Advance Wound Care Management and Infection Prevention with a vast array of specialized products.

The Sion Medical Group aim is to develop new innovative products that will influence, support and improve the quality of healthcare through the practice and management of infection control, wound management as well and the application of epidemiology in the healthcare settings.

The Company's Quality System is accredited to ISO 9001:2015, EN ISO 13485:2016, and in compliance with CFR Title 21 Parts 210, 211, 820. Our products are registered in the EU with CE certification and in the US under the FDA.



Povidone-Iodine

Povidone-Iodine products are antiseptic for topical application.

The product is used as a:

1. First aid antiseptic - Helps prevent skin infection in minor cuts, scrapes, and burns.
2. Preoperative patient skin preparation - Helps reduce bacteria that potentially can cause skin infection.

The products are packaged in convenient bottles, and easy-to-tear pouches.

Item No.	Description	Pkg.
54W202-1406-000	PVP Prep Pad	100 pc/bx, 10 bx/cs
54W402-4215-000	PVP Swabstick 1	50 pc/bx, 10 bx/cs
54W402-4234-000	PVP Swabsticks 3's	25 pc/bx, 10 bx/cs
54W413-4215-000	PVP Gel Swabstick 1	50 pc/bx, 10 bx/cs
307721	Prep Solution 2 oz	48 bottles/cs
307741	Prep Solution 3 oz	36 bottles/cs
307751	Prep Solution 4 oz	48 bottles/cs
307781	Prep Solution 8 oz	24 bottles/cs
307761	Prep Solution 16 oz	12 bottles/cs
307732	Prep Solution 32 oz	12 bottles/cs
307791	Prep Solution Gallon oz	4 bottles/cs
307821	Scrub Solution 2 oz	48 bottles/cs
307831	Scrub Solution 3 oz	36 bottles/cs
307841	Scrub Solution 4 oz	36 bottles/cs
307881	Scrub Solution 8 oz	24 bottles/cs
307861	Scrub Solution 16 oz	12 bottles/cs
307832	Scrub Solution 32 oz	12 bottles/cs
307891	Scrub Solution Gallon	4 bottles/cs



*Available in selected countries, for more information please contact info@sn-medical.com.

Alcohol

For topical application to clean the intended area before injection, venipuncture, or any form of medical procedure. Helps prevent the risk of skin infection in minor cuts, scrapes, and burns.

- The laminate provides hermetic closure of the product to prevent leakage and drying.
- Impregnated with 70% V/V Isopropyl Alcohol for optimal antibacterial action.
- Sterile prep pads are sterilized by Gamma irradiation.
- Latex-free.

Sterile

Item No.	Description	Pkg.
53W201-2208-000	Alcohol Prep Pad (L)	100 pc/bx, 10 bx/cs
53W201-1206-000	Alcohol Prep Pad (M)	100 pc/bx, 10 bx/cs



Non-Sterile

Item No.	Description	Pkg.
53W201-2206-000	Alcohol Prep Pad (L)	100 pc/bx, 10 bx/cs
54W201-1206-000	Alcohol Prep Pad (M)	100 pc/bx, 10 bx/cs
54W401-4215-000	Alcohol Swabsticks 1	50 pc/bx, 10 bx/cs
54W401-4234-000	Alcohol Swabsticks 3's	25 pc/bx, 10 bx/cs



Hydrogen Peroxide

Hydrogen Peroxide is intended for treatment on minor cuts and scrapes.

Item No.	Description	Pkg.
390061	Hydrogen Peroxide Bottle 16 oz	24 bottles/cs
307701	Hydrogen Peroxide Bottle 32 oz	12 bottles/cs

*Available in selected countries, for more information please contact info@sn-medical.com.

Protective Barrier Wipes

Protective Barrier Wipes is a liquid film-forming dressing that, upon application to intact skin, forms a protective film to help reduce friction during removal of tapes and films. It leaves a thin breathable film on the skin, which protects it against potential damage. It dries within seconds and does not build up residue on the skin, leaving it ready to apply with a new barrier.

The wipe's unique packaging makes it practical for traveling or use on-the-go. The product is according to FDA and CE (MDR) regulation.

Item No.	Description	Pkg.
54W211-1205-000	Protective Wipes	50 pc/bx, 40 bx/cs
54W211-1206-000	Protective Wipes	100 pc/bx, 20 bx/cs



Adhesive Remover Wipes

Adhesive Remover Wipes is a uniquely formulated solution which ensures pain-free removal of medical adhesive products. It is a non-invasive medical device that aids removal from intact skin, including:

1. Adhesive Tapes
2. Wound Dressings
3. Stoma Pouches

The product is according to FDA and CE (MDR) regulation.

Item No.	Description	Pkg.
54W209-1205-000	Adhesive Wipes	50 pc/bx, 40 bx/cs
54W209-1206-000	Adhesive Wipes	100 pc/bx, 20 bx/cs



*Available in selected countries, for more information please contact info@sn-medical.com.

Sting Free

Sting Free swabsticks and bottle are a liquid film-forming that upon application to intact skin, forms a protective film to help reduce friction during removal of tapes and films. It leaves a thin breathable film on the skin, which protects it against potential damage. It dries within seconds and does not build up residue on the skin, leaving it ready to apply with a new barrier.

The swabs's unique packaging makes it practical for traveling or use on-the-go.

Item No.	Description	Pkg.
54W411-4211-000	Swabsticks 1 ml	Contact Us
54W411-4213-000	Swabsticks 3 ml	Contact Us
	Spray Bottle 28 ml	Contact Us

Hand Sanitizer Wipes

Alcohol-based hand sanitizers effectively eliminate harmful viruses, fungi, and bacteria.

Item No.	Description	Pkg.
54W215-1000-000	Hand Sanitizer Wipes	1000 units/cs



Skin Electrode Pads

Disposable, easy-to-use, adhesive skin electrode pads.

Item No.	Description	Pkg.
54W210-1206-000	Skin Electrode Pads	100 pc/bx, 10 bx/cs



*Available in selected countries, for more information please contact info@sn-medical.com.

Smartphone Cleansing Wipe

- Travel friendly, cleans, removes oils and fingerprints.
- Individually wrapped wipe, protects smartphones from germs.
- Contains Isopropyl Alcohol USP, 40% v/v.ene.

Item No.	Description	Pkg.
54W203-6666-000	Cleansing Wipe 12.5x15 cm	100 pc/bx, 12 bx/cs

*Available only in the USA, for more information please contact info@sn-medical.com.



Nail Polish Remover Pads

Nail polish remover pads for easy and quick removal of nail polish.

Item No.	Description	Pkg.
54W208-1206-000	Nail Polish Remover Pads	100 pc/bx, 10 bx/cs



*Available in selected countries, for more information please contact info@sn-medical.com.

Sterile Lubricant Jelly

For lubrication prior insertion of catheters, endoscopes, and surgical instruments. Will not harm surgical instruments, rubber or plastic. Intended for medical purposes to lubricate body orifices and to facilitate entry of diagnostic or therapeutic devices.

- Gamma sterilized.
- Water soluble.
- Non irritating to skin, tissue and mucous membranes.
- Non greasy.
- Non staining.
- Non sticky.

Item No.	Description	Pkg.
53J505-4209-000	Tube 42g	12 pc/bx, 6 bx/cs
53J505-8209-000	Tube 82g	1 pc/bx, 48 bx/cs
53J505-1009-000	Tube 113g (4 oz)	12 pc/bx, 6 bx/cs
53J505-1009-011	Tube 113g (4 oz)	1 pc/bx, 56 bx/cs
53J805-3011-000	Sachets 3g	144 pc/bx, 12 bx/cs
53J805-5001-000	Sachets 5g	144 pc/bx, 12 bx/cs
53J805-5002-000	Sachets 5g	Bulk
53J805-6006-000	Sachets 15g	36 pc/bx, 144 bx/cs
53J505-1006-000	Syringes 6g	10 pc/bx, 120 bx/cs
53J505-1001-000	Syringes 8g	10 pc/bx, 120 bx/cs
53J505-1010-000	Syringes 10g	10 pc/bx, 120 bx/cs



Non-Sterile Lubricant Jelly

Item No.	Description	Pkg.
54J505-4209-000	Tube 42g	12 pc/bx, 6 bx/cs
54J505-8209-000	Tube 82g	1 pc/bx, 48 bx/cs
54J505-1009-000	Tube 113g (4 oz)	12 pc/bx, 6 bx/cs
54J505-1009-011	Tube 113g (4 oz)	1 pc/bx, 56 bx/cs
54J805-3011-000	Sachets 3g	144 pc/bx, 12 bx/cs
54J805-5001-000	Sachets 5g	144 pc/bx, 12 bx/cs
54J805-5002-000	Sachets 5g	Bulk
54J805-6006-000	Sachets 15g	36 pc/bx, 144 bx/cs

*Available in selected countries, for more information please contact info@sn-medical.com.

Cervilube™

Carbomer, Paraben & CHG Free.

Validated by Hologic for use with the ThinPrep® Pap test when used as instructed.

See www.hologic.com/thinprep for more details.

CerviLube™ Lubricating Jelly is intended for medical purposes, which is used to lubricate a body orifice to facilitate entry of a diagnostic or therapeutic device.

Advantages:

- Tested safe and effective for use with Pap Test.
- Eases insertion of diagnostic and therapeutic devices.
- Transparent and greaseless.
- Sterile/Non-sterile options.
- FDA and CE regulation.
- Complies with California Prop 65, unlike other lubricants commonly used by OB/GYN offices.

CerviLube™ Lubricating Jelly does not contain:

- Carbomers - Which impact the results of the PapTest.
- Parabens - Linked to increased risk of breast cancer.
- Chlorhexidine gluconate (CHG) - Possible skin irritation and allergic reactions.

Sterile

Item No.	Description	Pkg.
53J506-4209-000	Tube 42g	12 pc/bx, 6 bx/cs
53J506-8209-000	Tube 82g	1 pc/bx, 48 bx/cs
53J506-8200-000	Tube 82g (3 oz)	12 pc/bx, 6 bx/cs
53J506-1009-000	Tube 113g (4 oz)	12 pc/bx, 6 bx/cs
53J806-3012-000	Sachet 3g	144 pc/bx, 12 bx/cs
53J806-5011-000	Sachet 5g	144 pc/bx, 12 bx/cs



Non-Sterile

Item No.	Description	Pkg.
54J506-4209-000	Tube 42g	12 pc/bx, 6 bx/cs
54J506-8209-000	Tube 82g	1 pc/bx, 48 bx/cs
54J506-8200-000	Tube 82g (3 oz)	12 pc/bx, 6 bx/cs
54J506-1009-000	Tube 113g (4 oz)	12 pc/bx, 6 bx/cs
54J806-3012-000	Sachet 3g	144 pc/bx, 12 bx/cs
54J806-5011-000	Sachet 5g	144 pc/bx, 12 bx/cs



*Available in selected countries, for more information please contact info@sn-medical.com.

White Petrolatum

White Petrolatum is a skin protectant that soothes the skin by creating a protective moisture barrier.

- Contains white petrolatum, USP, 100% w/w.
- Does not contain artificial flavors, parabens, or other preservatives.
- Provides visible relief of dry skin when used regularly.
- Helps treat and prevent diaper rash.
- Temporarily protects minor cuts, scrapes, and burns.
- Temporarily helps protect and relieve chapped and cracked skin or lips.

Item No	Description	Pkg.
54J510-0010-000	Tube 10g (0.35 oz)	50 pc/bx, 12 bx/cs
54J810-5201-000	Sachet 5g (0.18 oz)	144 pc/bx, 6 bx/cs



Vitamins A&D Ointment

Vitamins A&D works by moisturizing and sealing the skin. Ideal for the treatment of any of the following conditions: diaper rash, chapped and cracked skin or lips, minor cuts and burns.

Item No	Description	Pkg.
54J512-0419-011	Tube 113g (4 oz)	1 pc/bx, 56 bx/cs
54J512-0419-000	Tube 113g (4 oz)	12 pc/bx, 6 bx/cs
54J811-5201-000	Sachet 5g (0.18 oz)	144 pc/bx, 6 bx/cs



*Available in selected countries, for more information please contact info@sn-medical.com.

Antibiotics Products

The antibiotic products are divided into two different products: Bacitracin Zinc Ointment and Triple Antibiotic Ointment. This antibiotic ointment is a must-have, over-the-counter product, in your first aid kit. It is used for topical treatment of a variety of localized skin infections in minor scrapes, cuts and burns, as well as for the prevention of wound infections.

Bacitracin Zinc Ointment

Item No	Description	Pkg.
54C810-0509-000	Tube 14g (0.5 oz)	100 pc/cs - Bulk
54C810-0109-000	Tube 28g (1 oz)	200 pc/cs - Bulk
54C810-9201-000	Sachet 10g (0.03 oz)	144 pc/bx, 12 bx/cs



Triple Antibiotic Ointment

Item No	Description	Pkg.
54C812-0509-000	Tube 14g (0.5 oz)	100 pc/cs - Bulk
54C812-0109-000	Tube 28g (1 oz)	200 pc/cs - Bulk
54C812-9201-000	Sachet 10g (0.03 oz)	144 pc/bx, 12 bx/cs



*Available in selected countries, for more information please contact info@sn-medical.com.

Hydrocortisone

Hydrocortisone Acetate Cream is an over-the-counter drug. It is used for the temporary relief of itching associated with minor skin irritations: allergy, inflammation, rashes due to eczema, insect bites, poison ivy, poison oak, poison sumac, detergents, cosmetics, jewelry, seborrheic dermatitis, psoriasis, external genital, feminine, and anal itching.

Item No	Description	Pkg.
54C813-0509-000	Tube 14g (0.5 oz)	100 pc/cs - Bulk
54C813-0109-000	Tube 28g (1 oz)	200 pc/cs - Bulk
54C813-9401-000	Sachet 10g (0.03 oz)	144 pc/bx, 12 bx/cs



*Available in selected countries, for more information please contact info@sn-medical.com.

Unna Boot Flexible Bandages

With Zinc Paste & Calamine.

The Unna Boot Bandage is ideal for the treatment and management of any of the following diagnosed conditions:

- Venous insufficiency.
- Lower extremity edema.
- Venous stasis ulceration.
- Lymphedema.
- Minor fractures, sprains and strains.
- Skin grafts.

Unna Boot is available with and without Calamine.

Item No.	Description	Pkg.
54C257-7590-000	Zinc Paste & Zinc Calamine 3"	1 pc/bx, 12 bx/cs
54C257-1090-000	Zinc Paste & Zinc Calamine 4"	1 pc/bx, 12 bx/cs
54C156-7590-000	Zinc Paste 3"	1 pc/bx, 12 bx/cs
54C156-1090-000	Zinc Paste 4"	1 pc/bx, 12 bx/cs



Parassan®

- Non-woven polyester and viscose pads indicated for superficial wounds such as contusions, abrasions, as well as various other dermatological conditions.
- Keeping skin wounds moist.
- Barrier between the wound and the surrounding environment.
- Prevents adhesion between the wound and the dressing.

Parassanet Pads

- Parassanet is the perfect replacement of impregnated gauzes and is manufactured from 100% cotton, this is ideal for the treatment of burn wounds.
- The open weave net allows for the passage of wound exudates to be absorbed by a secondary dressing, avoiding maceration of the surrounding tissue.
- Non adhering to the wound.
- Provides a moist healing environment.

*Available in selected countries, for more information please contact info@sn-medical.com.

Item No.	Description	Pkg.
105012	Parassan® Pad 10x10 cm (4"x4")	10 pc/bx, 24 bx/cs
106012	Parassan® Pad 10x20 cm (4"x8")	10 pc/bx, 12 bx/cs
105112	Parassanet Pad 10x10 cm (4"x4")	10 pc/bx, 24 bx/cs
106112	Parassanet Pad 10x20 cm (4"x8")	10 pc/bx, 12 bx/cs
53C600-1090-000	Parassan® Bandage 10x7 m (4" x 7.66yd)	1 pc/bx, 12 bx/cs



Medical Grade Honey

L-Mesitran Medical-Grade Honey products cure complex wounds, remove biofilms from infected wounds, and prevent infections and antibiotic-resistance problems, making healing faster and more efficient.

L-Mesitran Soft | L-Mesitran Ointment

L-Mesitran Soft and Ointment contain medical-grade honey, suitable for the treatment of burns, chronic, oncological, and acute wounds such as abrasions and cuts.

L-Mesitran Soft

Item No.	Description	Pkg.
54J515-5015-130	Tube 15 g	1 pc/bx, 14 bx/cs
54J515-4050-101	Tube 50 g	1 pc/bx, 12 bx/cs
53J909-0804-000	Syringe 8 g Sterile	4 pc/bx, 70 bx/cs
53J909-1201-000	Syringe 5 g Sterile	12 pc/bx, 24 bx/cs
53J909-1001-000	Syringe 10 g Sterile	12 pc/bx, 20 bx/cs
54J615-3250-101	Jar 250 g	30 pc/cs
54J808-5009-726	Sachet 0.9 g	10 pc/bx, 72 bx/cs

L-Mesitran Ointment

Item No.	Description	Pkg.
54J715-0020-726	Tube 20 g	1 pc/bx, 14 bx/cs
54J715-0050-726	Tube 50 g	1 pc/bx, 12 bx/cs
54J715-0250-726	Jar 250 g	1 pc/bx, 30 bx/cs



*Available in selected countries, for more information please contact info@sn-medical.com.

L-Mesitran Foam | Tulle | Border | Net | Hydro | Active

These low-adherence wound dressing products contain medical-grade honey. They are suitable for treating and protecting a wide range of wounds such as superficial, chronic, acute wounds, cuts, surgical wounds, and burns.

Item No.	Description	Pkg.
53J904-1108-726	Border 10x10 cm (Pad)	10 pc/bx, 20 bx/cs
53J904-1107-726	Active 7.5x4.5 cm (Pad)	5 pc/bx, 72 bx/cs
53J904-1110-726	Hydro 10x10 cm (Pad)	10 pc/bx, 20 bx/cs
109116	Foam	10 pc/bx, 20 bx/cs
108218	Net (Pouch)	10 pc/bx, 20 bx/cs
108124	Tulle (Pouch)	10 pc/bx, 20 bx/cs



Hydrogels

Amorphous Hydrogel Tubes

The Amorphous Hydrogel is Carbomer, Paraben & CHG free. Not made with natural rubber latex, non-cytotoxic, non-irritating. Multipurpose, non-sterile product, intended to maintain a moist area by providing moisture to dry or minimally exudating wounds and injuries.

Recommended Use:

- Acute wounds.
- Burns.
- First- and second-degree burns.
- Leg ulcers.
- Non/minimally exudating wounds.
- Partial - and full-thickness wounds.
- Pressure ulcers.
- Superficial wounds.
- Surgical wounds.
- Venous ulcers.
- Lacerations, abrasions, and skin tears.

Item No.	Description	Pkg.
54J506-0109-000	Tube 1 oz - Nuzzle	1 pc/bx
54J506-0309-000	Tube 3 oz - Flip top	1 pc/bx



*Available in selected countries, for more information please contact info@sn-medical.com.

Sion HS-Spirit

Post Operative Antibacterial Hydrogel Dressing for the purpose of treating surgical site infection. HS-Spirit is an antibacterial barrier, intended to create an ideal wound environment and stimulates wound healing for acute and chronic wounds. Transparent, self-adhesive borders, waterproof, and does not adhere to the wound.



Sion MIS HS-Spirit

Post Operative Antibacterial Hydrogel Dressing for minimal invasive surgery. Antibacterial barrier, for (superficial) acute and chronic wounds. Creates an ideal wound environment and stimulates wound healing. Transparent, self-adhesive borders (waterproof), does not adhere to the wound. Can stay on the wound for 3-5 days.

Item No.	Description	Pkg.
53J904-1037-000	Hydrogel Dressing (HS-Spirit) 10x37 cm	10 pc/bx, 6 bx/cs
53J904-1037-011	Hydrogel Dressing (HS-Spirit) 10x37 cm	1 pc/bx
53J904-1107-000	Hydrogel Dressing (MIS) 4.5x7.5 cm	5 pc/bx, 12 bx/cs



*Available in selected countries, for more information please contact info@sn-medical.com.

Gauze Pads

Gauze Pads made of 100% cotton, for cleaning, dressing, prepping wounds. Offered as sterile, non-sterile, x-ray gauze pads, and Y-cut gauze pads.

Sterile Gauze Pads

Item No	Description	Pkg.
001012	Gauze Pads 5x5 cm 8 ply 17 Thread	1 pc/pck, 100 pck/bx, 52 bx/cs
001042	Gauze Pads 5x5 cm 8 ply 17 Thread	4 pcs/pck, 375 pck/cs
010012	Gauze Pads 7.5x7.5 cm 8 ply 17 Thread	1 pc/pck, 100 pck/bx, 30 bx/cs
010022	Gauze Pads 7.5x7.5 cm 8 ply 17 Thread	2 pcs/pck, 750 pck/cs
010032	Gauze Pads 7.5x7.5 cm 8 ply 17 Thread	2 pcs/pck, 600 pck/cs
010042	Gauze Pads 7.5x7.5 cm 8 ply 17 Thread	4 pcs/pck, 200 pck/cs
010052	Gauze Pads 7.5x7.5 cm 8 ply 17 Thread	5 pcs/pck, 150 pck/cs
010092	Gauze Pads 7.5x7.5 cm 8 ply 17 Thread	10 pcs/pck, 40 pck/cs
013052	Gauze Pads 7.5x7.5 cm 12 ply 17 Thread	5 pcs/pck, 130 pck/cs
030012	Gauze Pads 10x10 cm 8 ply 17 Thread	1 pc/pck, 100 pck/bx, 20 bx/cs
030022	Gauze Pads 10x10 cm 8 ply 17 Thread	2 pcs/pck, 500 pck/cs
030052	Gauze Pads 10x10 cm 8 ply 17 Thread	5 pcs/pck, 72 pck/cs
032022	Gauze Pads 10x10 cm 12 ply 17 Thread	2 pcs/pck, 400 pck/cs
032052	Gauze Pads 10x10 cm 12 ply 17 Thread	5 pcs/pck, 60 pck/cs
032092	Gauze Pads 10x10 cm 12 ply 17 Thread	10 pcs/pck, 15 pck/cs
032292	Gauze Pads 10x10 cm 12 ply 17 Thread Double Pack	10 pcs/pck, 10 pck/cs
034292	Gauze Pads 10x10 cm 16 ply 17 Thread Double Pack	10 pcs/pck, 10 pck/cs
300012	Gauze Eye Pad 1pc Sterile	1 pc/pck, 70 pck/bx, 24 bx/cs

*Available globally, for more information please contact info@sn-medical.com.

Non-Sterile Gauze Pads in Bulk

Item No	Description	Pkg.
009000	Gauze Pads 7.5x7.5 cm 8 ply 13 Thread	100 pcs/bag, 50 bags/cs
010000	Gauze Pads 7.5x7.5 cm 8 ply 17 Thread	100 pcs/bag, 50 bags/cs
013000	Gauze Pads 7.5x7.5 cm 12 ply 17 Thread	100 pcs/bag, 40 bags/cs
029000	Gauze Pads 10x10 cm 8 ply 13 Thread	100 pcs/bag, 30 bags/cs
030000	Gauze Pads 10x10 cm 8 ply 17 Thread	100 pcs/bag, 25 bags/cs
032000	Gauze Pads 10x10 cm 12 ply 17 Thread	100 pcs/bag, 20 bags/cs
034000	Gauze Pads 10x10 cm 16 ply 17 Thread	100 pcs/bag, 20 bags/cs
044002	Gauze Pads 10x20 cm 16 ply 17 Thread	100 pcs/bag, 8 bags/cs



Gauze Pads

X-Ray Gauze Pads

Item No	Description	Pkg.
013294	Gauze Pads 7.5x7.5 cm 12 ply 17 Thread with X-Ray Sterile Double Pack	10 pcs/pck, 200 pck/cs
032294	Gauze Pads 10x10 cm 12 ply 17 Thread with X-Ray Sterile Double Pack	10 pcs/pck, 100 pck/cs
034294	Gauze Pads 10x10 cm 16 ply 17 Thread with X-Ray Sterile Double Pack	10 pcs/pck, 100 pck/cs
013001	Gauze Pads 7.5x7.5 cm 12 ply 17 Thread with X-Ray Non-Sterile	100 pcs/bag, 40 bags/cs
032001	Gauze Pads 10x10 cm 12 ply 17 Thread with X-Ray Non-Sterile	100 pcs/bag, 20 bags/cs

Y-Cut Gauze Pads

Item No	Description	Pkg.
010612	Gauze Pads with Y-Cut 7.5x7.5 cm 8 ply Sterile	1 pc/pck, 100 pck/bx, 24 bx/cs
010622	Gauze Pads with Y-Cut 7.5x7.5 cm 8 ply Sterile	2 pcs/pck, 1500 pck/cs

*Available globally, for more information please contact info@sn-medical.com.

Gauze Balls

Highly absorbent Gauze Balls for multiple uses, offered in a variety of sizes and categories.

Item No	Description	Pkg.
063032	Gauze Balls type G-3 (2 ply) 35mm Sterile	3 pcs/pck, 360 pck/cs
063042	Gauze Balls type G-3 (2 ply) 35mm Sterile	4 pcs/pck, 225 pck/cs
067052	Gauze Balls type H-4 (1 ply) 45mm Sterile	5 pcs/pck, 100 pck/cs
063062	Gauze Balls type G-3 (2 ply) 35mm Sterile	6 pcs/pck, 500 pck/cs
063002	Gauze Balls type G-3 (2 ply) 35mm Bulk	500 pcs/bag, 8 bags/cs
064002	Gauze Balls type G-4 (2 ply) 45mm Bulk	500 pcs/bag, 6 bags/cs
069000	Gauze Balls type K-4 (2 ply) 45mm Bulk	500 pcs/bag, 5 bags/cs
067000	Gauze Balls type H-4 (1 ply) 45mm Bulk	500 pcs/bag, 10 bags/cs
064001	Gauze Balls G-4 XR (2 ply) 45mm Bulk	500 pcs/bag, 6 bags/cs
110294	Peanuts Balls 8x8 cm with XR Sterile Double Pack	10 pcs/pck, 30 pck/cs
750702	Cotton Balls 0.5g 100 pcs Bag Sterile	100 pcs/bag, 50 bags/cs



Gauze Balls



Gauze Balls XR

Non-Woven Pads

Non-Woven Pads for fluid and blood absorption. Does not leave any fiber or residue on the body.

Item No	Description	Pkg.
090000	Non-Woven Pads 5x5-4 ply V70/P30 30gsm	100 pcs/bag, 50 bags/cs
090600	Non-Woven Pads 5x5-6 ply V70/P30 30gsm	100 pcs/bag, 50 bags/cs
092000	Non-Woven Pads 7.5x7.5-4 ply V70/P30 30gsm	100 pcs/bag, 40 bags/cs
092600	Non-Woven Pads 7.5x7.5-6 ply V70/P30 30gsm	100 pcs/bag, 40 bags/cs
094000	Non-Woven Pads 10x10-4 ply V70/P30 30gsm	100 pcs/bag, 32 bags/cs
094600	Non-Woven Pads 10x10-6 ply V70/P30 30gsm	100 pcs/bag, 32 bags/cs



Non-Woven Pads 7.5x7.5 - 4 Ply

*Available globally, for more information please contact info@sn-medical.com.

Non-Woven Balls

Non-Woven Balls are highly absorbent and are used for skin prepping, wound cleansing, and for the use of general-purpose.

Item No	Description	Pkg.
072000	Non-Woven Balls H2 20mm V70/P30 30gsm	500 pcs/bag, 10 bags/cs
073000	Non-Woven Balls H3 30mm V70/P30 30gsm	500 pcs/bag, 8 bags/cs
074000	Non-Woven Balls H4 40mm V70/P30 30gsm	500 pcs/bag, 6 bags/cs

Bandages

Gauze Bandage

Offering a variety of Bandages for the use of wound and limb wrapping.

Item No	Description	Pkg.
051805	Gauze Bandage 2.7m x 5cm	10 pcs/pck, 84 pck/cs
051807	Gauze Bandage 2.7m x 7cm	10 pcs/pck, 84 pck/cs
051810	Gauze Bandage 2.7m x 10cm	10 pcs/pck, 40 pck/cs
051812	Gauze Bandage 2.7m x 12cm	10 pcs/pck, 40 pck/cs



Gauze Bandage 2.7m x 10cm

Conforming Stretch Bandage

Item No	Description	Pkg.
54B111-4525-100	Conforming Bandage 1000 4.5cm x 2.7m	5 pcs/pck, 20 pck/cs
54B110-5025-000	Conforming Bandage 1000 5cm x 2.7m	5 pcs/pck, 40 pck/cs
059205	Conforming Bandage 1000 5cm x 2.7m	5 pcs/pck, 200 pck/cs
54B111-6525-100	Conforming Bandage 1000 6.5cm x 2.7m	5 pcs/pck, 20 pck/cs
059208	Conforming Bandage 1000 7.5cm x 2.7m	5 pcs/pck, 200 pck/cs
54B110-7525-000	Conforming Bandage 1000 7.5cm x 2.7m	5 pcs/pck, 40 pck/cs
54B110-9025-100	Conforming Bandage 1000 9cm x 2.7m	5 pcs/pck, 20 pck/cs
059210	Conforming Bandage 1000 10cm x 2.7m	5 pcs/pck, 160 pck/cs
54B110-1025-000	Conforming Bandage 1000 10cm x 2.7m	5 pcs/pck, 20 pck/cs
54B110-1125-100	Conforming Bandage 1000 11cm x 2.7m	5 pcs/pck, 20 pck/cs



Conforming Bandage

*Available globally, for more information please contact info@sn-medical.com.

Elastic Bandage

Item No	Description	Pkg.
54B661-6036-100	Elastic Bandage 600 6cm x 3.7m	1 pc, 360 pcs/cs
54B660-6041-000	Elastic Bandage 600 6cm x 4.5m + Clip	1 pc, 250 pcs/cs
54B661-8036-100	Elastic Bandage 600 8cm x 3.7m	1 pc, 280 pcs/cs
54B660-8040-000	Elastic Bandage 600 8cm x 4.57m	1 pc, 250 pcs/cs
54B660-8041-000	Elastic Bandage 600 8cm x 4.5m + Clip	1 pc, 250 pcs/cs
54B661-1036-100	Elastic Bandage 600 10cm x 3.7m	1 pc, 240 pcs/cs
54B660-1040-000	Elastic Bandage 600 10cm x 4.57m	1 pc, 250 pcs/cs
54B660-1041-000	Elastic Bandage 600 10cm x 4.5m + Clip	1 pc, 250 pcs/cs
54B661-1236-100	Elastic Bandage 600 12cm x 3.7m	1 pc, 200 pcs/cs
54B661-1536-100	Elastic Bandage 600 15cm x 3.7m	1 pc, 160 pcs/cs
54B660-1241-000	Elastic Bandage 600 12cm x 4.5m + Clip	1 pc, 250 pcs/cs
54B660-1541-000	Elastic Bandage 600 15cm x 4.5m + Clip	1 pc, 100 pcs/cs
54B660-1540-000	Elastic Bandage 600 15cm x 4.57m	1 pc, 100 pcs/cs
54B675-8040-000	Elastic Bandage 750 8cm x 4.5m	1 pc, 250 pcs/cs
54B675-8041-000	Elastic Bandage 750 8cm x 4.57m + Clip	1 pc, 250 pcs/cs
54B675-1040-000	Elastic Bandage 750 10cm x 4.5m	1 pc, 250 pcs/cs
54B675-1041-000	Elastic Bandage 750 10cm x 4.57m + Clip	1 pc, 250 pcs/cs
54B675-1240-000	Elastic Bandage 750 12cm x 4.5m	1 pc, 250 pcs/cs
54B675-1273-100	Elastic Bandage 750 12cm x 7m + Clip	1 pc, 100 pcs/cs
54B675-1540-000	Elastic Bandage 750 15cm x 4.5m	1 pc, 100 pcs/cs
54B675-1573-100	Elastic Bandage 750 15cm x 7m + Clip	1 pc, 100 pcs/cs



Elastic Bandage

Fluffy Gauze

Absorbent Advanced Bandages for wound and limb wrapping to protect injuries.

Item No	Description	Pkg.
53B990-1130-000	Multiply Fluffy Gauze Bandage 11.5cm x 3.75m Sterile	1 pc/pck, 100 pck/cs
53B990-1133-000	Multiply Fluffy Gauze Bandage 11.5cm x 3.75m XR Sterile	1 pc/pck, 100 pck/cs



Multiply Bandage

*Available globally, for more information please contact info@sn-medical.com.

Laparotomy Sponges

Laparotomy Sponges for the use during operations and surgical procedures, keeping the surgical area free of excessive body fluids, bloods and other liquids.

Sterile

Item No	Description	Pkg.
231254	Laparotomy Sponges 30x30 cm 4 ply Loop + Washed with X-Ray Double Pack	5 pcs/pck, 60 pck/cs
233254	Laparotomy Sponges 30x30 cm 6 ply Loop + Washed with X-Ray Double Pack	5 pcs/pck, 60 pck/cs
234254	Laparotomy Sponges 35x35 cm 6 ply Washed with X-Ray Double Pack	5 pcs/pck, 50 pck/cs
239254	Laparotomy Sponges 30x30 cm 8 ply Loop + Washed with X-Ray Double Pack	5 pcs/pck, 50 pck/cs
261254	Laparotomy Sponges 45x45 cm 4 ply Loop with X-Ray Double Pack	5 pcs/pck, 50 pck/cs
271224	Laparotomy Sponges 45x45 cm 4 ply Washed with X-Ray Double Pack	2 pcs/pck, 120 pck/cs
271254	Laparotomy Sponges 45x45 cm 4 ply Washed with X-Ray Double Pack	5 pcs/pck, 50 pck/cs
281254	Laparotomy Sponges 45x45 cm 4 ply Loop + Washed with X-Ray Double Pack	5 pcs/pck, 50 pck/cs
289254	Laparotomy Sponges 45x45 cm 8 ply Loop + Washed with X-Ray Double Pack	5 pcs/pck, 20 pck/cs

Non-Sterile

Item No	Description	Pkg.
231000	Laparotomy Sponges 30x30 cm 4 ply Loop + Washed with X-Ray	600 pcs/cs
239001	Laparotomy Sponges 30x30 cm 8 ply Loop + Washed with X-Ray	400 pcs/cs
281000	Laparotomy Sponges 45x45 cm 4 ply Loop + Washed with X-Ray	400 pcs/cs
289000	Laparotomy Sponges 45x45 cm 8 ply Loop + Washed with X-Ray	200 pcs/cs

*Available globally, for more information please contact info@sn-medical.com.

Surgical Towels (OR Towels)

Disposable Surgical Towels made of 100% de-linted cotton, used for patient infection prevention, absorption of small amounts of fluids and general use in the operation room environment. Available with or without sterilization paper wrap. Sterile and non-sterile presentations. With and without detectable X-ray thread.

Main features:

- Very low linting.
- Pre-washed.
- Highly absorbent.
- Corner-fold presentation for easy use.
- Latex free.

Description	Pkg.
Surgical Towel 40x66 cm (16"x26") Blue Sterile	1 pcs/pck
Surgical Towel 40x66 cm (16"x26") Blue Sterile	2 pcs/pck
Surgical Towel 40x66 cm (16"x26") Blue Sterile	4 pcs/pck
Surgical Towel 40x66 cm (16"x26") Blue Sterile	6 pcs/pck
Surgical Towel 40x66 cm (16"x26") Blue Sterile	8 pcs/pck
Surgical Towel 40x66 cm (16"x26") Blue Sterile	10 pcs/pck
Surgical Towel 40x66 cm (16"x26") Blue Non-Sterile	100 in Bulk
Surgical Towel 40x66 cm (16"x26") Blue Non-Sterile	200 in Bulk
Surgical Towel 40x66 cm (16"x26") Blue Non-Sterile	400 in Bulk



*Available globally, for more information please contact info@sn-medical.com.

Absorbent Pad

Non-Woven Absorbent Pad


Item No	Description	Pkg.
144012	Non-Woven ABD Pad 10x10 cm Sterile	1 pc/pck, 500 pck/cs
144212	Non-Woven ABD Pad 10x20 cm Sterile	1 pc/pck, 400 pck/cs
144612	Non-Woven ABD Pad 20x20 cm Sterile	1 pc/pck, 300 pck/cs
144712	Non-Woven ABD Pad 20x25 cm Sterile	1 pc/pck, 250 pck/cs
144812	Non-Woven ABD Pad 20x30 cm Sterile	1 pc/pck, 150 pck/cs
144912	Non-Woven ABD Pad 20x40 cm Sterile	1 pc/pck, 120 pck/cs
144000	Non-Woven ABD Pad 10x10 cm Non-Sterile	1200 pcs/cs
144200	Non-Woven ABD Pad 10x20 cm Non-Sterile	800 pcs/cs
144600	Non-Woven ABD Pad 20x20 cm Non-Sterile	600 pcs/cs
144700	Non-Woven ABD Pad 20x25 cm Non-Sterile	400 pcs/cs
144800	Non-Woven ABD Pad 20x30 cm Non-Sterile	300 pcs/cs
144900	Non-Woven ABD Pad 20x40 cm Non-Sterile	200 pcs/cs

Cotton Gauze Absorbent Pad

Item No	Description	Pkg.
134012	Gauze ABD Pad 10x10 cm Sterile	1 pc/pck, 500 pck/cs
134212	Gauze ABD Pad 10x20 cm Sterile	1 pc/pck, 400 pck/cs
134412	Gauze ABD Pad 10x35 cm Sterile	1 pc/pck, 300 pck/cs
134612	Gauze ABD Pad 20x20 cm Sterile	1 pc/pck, 250 pck/cs
134812	Gauze ABD Pad 20x30 cm Sterile	1 pc/pck, 150 pck/cs
134912	Gauze ABD Pad 20x40 cm Sterile	1 pc/pck, 120 pck/cs
134000	Gauze ABD Pad 10x10 cm Non-Sterile	1000 pcs/cs
134200	Gauze ABD Pad 10x20 cm Non-Sterile	800 pcs/cs
134400	Gauze ABD Pad 10x35 cm Non-Sterile	600 pcs/cs
134600	Gauze ABD Pad 20x20 cm Non-Sterile	400 pcs/cs
134800	Gauze ABD Pad 20x30 cm Non-Sterile	300 pcs/cs
134900	Gauze ABD Pad 20x40 cm Non-Sterile	200 pcs/cs

*Available globally, for more information please contact info@sn-medical.com.

Each set is available in different configurations and dimensions per the customer's request. Examples of sets appear on the following pages. Contact us for further specifications.



*Available globally, for more information please contact info@sn-medical.com.

Ophthalmic

Ophthalmic Pack

The ophthalmic set is suitable for eye surgeries such as cataract surgeries and pterygium surgeries. This set offers a variety of essential single-use items needed during most surgeries, such as pads and gauze, syringes and bowls, eye shields, and sort covers that are essential as a barrier between the sterile operation area and the contaminated area.

Description of component	Quantity
Gauze Ball G4	10
Sponge Forceps 19 cm	1
Tray 18x23 cm	1
Table Drape 75x100 cm	1
Applicator 6"	2
Cup 60ml	2
Polyethylene Bag 16x30 cm	1
Syringe 1LS	1
Needle 18GA	1
Syringe 10LL	1
Syringe 3LL	2
Syringe 5LL	1
Gauze Pad 8ply 7.5x7.5 cm	5
Section Irrigation Bag 30x50 cm	1
Eye Shield	1
Tray 13x17 cm	1
Eye Pad	1
Ophthalmic Drape 27.5x40 cm	1
Velcro Tube Holder 2.5x11.5 cm	1
Cup 60ml	2
Eye Applicator with Rubber	5
Cellulose Surgical Eye Spears	5
Image Cover 80x80 cm	2
Ophthalmic Drape 150x250 cm	1
Table Cover 100x140 cm	1
Standard Gown XL	2
Infusion Set Vented LS	1

*Available globally, for more information please contact info@sn-medical.com.

Orthopedic

Total Knee Replacement Pack

The set is intended for use in knee replacement surgeries. A complete set of covers is included in this set for all the equipment used in knee joint surgeries, such as a table cover and mayo cover, as well as patient covers to maintain a sterile environment during the operation, even when large amounts of fluid is involved.

Description of component	Quantity
Drape 180x200 cm with Adhesive	1
Mayo Cover 80x140 cm	2
Table Cover 137x220 cm	2
Velcro Tube Holder 2.5x11.5 cm	1
Foot Drape 170x200 cm	2
Stockinette S 30x75 cm	1
Split Sheet	1
Extremity Drape	1
O.P Tape 10x50 cm	1
Light Handle Cover	1
Reinf.Gown XL (SMS) + Towel	2

*Available globally, for more information please contact info@sn-medical.com.

General Surgery

Laparoscopy Pack

The Laparoscopy set is offered for minimally invasive pelvic cavity laparoscopy surgery. It offers a wide variety of medical equipment that is necessary for most procedures; including diathermia pencil, lap sponges, blades, camera cover, and other single-use covers such as gowns, table covers, mayo covers, and patient universal cover set that assures sterile operation area.

Description of component	Quantity
Lap Sponge 45x45 cm	10
Reinf. Gown XL + Towel	1
Reinf. Gown XXL + Towel	2
Suction Irrigation Bag 30x50 cm	1
Hand Towel 30x40 cm	2
Mayo Cover 60x190 cm	1
Table Cover 100x140 cm	1
Table Cover 152x250 cm	1
Camera Cover	1
Diathermy Pouch 2 Pocket 35x42 cm	1
Side Drape 75x100 cm	2
Foot Drape 170x200 cm	1
Head Drape 150x240 cm	1
Blade #11	1
Needle Counter	1
Diathermy Pencil 5 m	1
Suction Tube 350 cm	1
MIS HS-Spirit Medical Honey Dressing	1

*Available globally, for more information please contact info@sn-medical.com.

Cardiology

Cath Lab Pack

The Cardiovascular set is used for catheterization procedures. This set contains single-use products that are needed in every procedure, such as needles for the Selinger technique, Guidewire, syringes, needles, cups, blades, and a waste bowl used as a waste blood container. This set also includes covers that are essential as a barrier between the sterile operation area and the contaminated area.

Description of component	Quantity
Gauze Pad 12ply 10x10 cm	20
18GAX2.75 Introducer Needle	1
Guidewire 150 Length, 035 Diameter	1
Needle 18GA	1
Needle 21GA	1
Needle 25GA	1
Mini Scalpel #11	1
Blue Towel	4
Image Cover 80x80 cm	1
Protective Cover 80x100 cm	1
Radial Drape 280x350 cm	1
Table Cover 137x180 cm Without Reinf. Area	2
Syringe 3LL	2
Cup 60ml-Red	1
Cup 60ml-Green	1
Cup 60ml	2
Cup 100ml	2
Basin Guidewire 2500cc	1
Transparent Basin 250cc	1
Syringe 20LS	1
Syringe 20LL	1
Syringe 10LL	3
Waste Bowl	1
Sharps Holder	1
Blue Tray	1
Labels Set	1

*Available globally, for more information please contact info@sn-medical.com.

Angio Pack

The Angio set is used for peripheral angiography catheterization procedures. This set offers single-use products needed for most Angio procedures, including covers for keeping a sterile procedure environment, and other single-use products as in the next list.

Description of component	Quantity
Gauze Pad 12 ply 10x10 cm	10
Gauze Ball G3	10
Sponge Forceps 19 cm	1
Infusion Set Vented LL 190 cm	1
Needle 25GA	1
Stopcock 4 Way	3
Blue Tower For Surgical Room 40x66 cm	4
Image Cover 80x80 cm	2
Key Board Cover 40x60 cm	1
Protective Cover 80x100 cm	1
Angio Drape 240x350 cm	1
Table Cover 137x180 cm	1
Needle 18GA	1
Needle 23G	1
Syringe 3LL	1
Scalpel #15	1
Needle Holder	1
Scalpel #11	1
Kocher Forceps	4
Extension 50 cm	1
Basin Guidewire 2500cc	2
Transparent Basin 250cc	2
Blue Basin 500cc	1
Syringe 20LL	3
Syringe 10LL	4
Sharps Holder	1
Kidney Basin 700cc	1

*Available globally, for more information please contact info@sn-medical.com.

Neurosurgery Set

The Neurosurgery set is used in neurosurgery procedures. The set offers all covers needed for tables and other equipment needed for this operation. The patient drape has a transparent area for a better view of the surrounding area and a wide fluid pouch around the patient's head. This set contains sponges, laps, and other essential unique products such as suction tubes and skin staplers.

Description of component	Quantity
Gauze Pad 12 ply 10x10 cm XR	20
Lap Sponge 45x45 cm	10
Reinf. Gown XL + Towel	1
Skin Marker - Standard Tip	4
Tip Cleaner 4.5x4.5 cm	1
Light Handle Cover	6
Tray 19.5x22 cm	1
Blue Towel For Surgical Room 40x66 cm	8
Neuro Transparent Drape 220x240 cm	1
Mayo Cover 60x190 cm	1
Table Cover 100x140 cm	2
Table Cover 152x250 cm	1
Diathermy Pouch 2 Pocket 35x42 cm	2
Table Drape 75x100 cm	1
Line 14.5 cm	2
Needle 23g	2
Tommy Syringe 50cc	2
Syringe 10 LL	2
Syringe 20 LS	2
Cup 100cc	4
Scalpel #11	1
Blade #10	2
Skin Stapler	1
Needle Counter	1
Diathermy Pencil 5 m	1
Suction Tube 350 cm	3
Wide Sponge 6x6 cm	9
Venflon 20g	3

*Available globally, for more information please contact info@sn-medical.com.

Lemon - Glycerin Swabsticks

Lemon Glycerin Swabsticks for oral care are ideal for those needing extra moistening due to medication side effects, radiation therapy, dental surgery, disease, after general anesthesia, and drug therapy. In addition, they are great for patients prescribed fluid restricted diets. They stimulate the saliva glands, providing refreshment and relief for patients with restricted oral intake.

- Pre-moistened with refreshing lemon-flavored glycerine-based dosage form.
- Sugar-free.
- For temporary relief of minor discomfort and protection of irritated areas in soremouth and sore throat.
- Disposable. Contain three single-use swabsticks in each foil pack.

Item No	Description	Pkg.
54W403-4234-000	10cm (4") 3 Swabs per Sachet	25 pc/bx, 10 bx/cs
54W403-6234-000	15cm (6") 3 Swabs per Sachet	25 pc/bx, 10 bx/cs



Ice Cool Bandage

Ice Cool is an elastic cold bandage that combines cooling and compression. It offers instant pain relief and cools for hours bruises, inflammations, and sprains.

Item No	Description	Pkg.
54C668-8026-000	Elastic Cool Bandage 8cm x 4m	1 pc/bx, 12 bx/cs



*Available globally, for more information please contact info@sn-medical.com.

Nissoft Towels

Multi-purpose Nissoft Towels are made of non-woven fabric and provided in convenient packaging for retrieval.

Item No	Description	Pkg.
098200	Towel 20x20 cm	200 pcs/bx, 12 bx/cs
098300	Towel 30x30 cm	24 pcs/bx, 20 bx/cs
098400	Towel 30x40 cm	50 pcs/bx, 18 bx/cs
098500	Towel 40x40 cm	40 pcs/bx, 12 bx/cs



Baby Blankets

Multipurpose Baby Blankets for indoor and outdoor use, perfect as a burping cloth and very convenient as a changing mat anywhere. Great as a sheet cover for the Moses Basket and Cot.

- Single use - No washing required.
- Ultra absorbent, hypoallergenic, original unprinted.
- Heat retaining.
- Very smooth and comfortable.

Color: White (additional colors available per MOQ).

Item No	Description	Pkg.
098212	72x60 cm - Sterile	10 pcs, 25 pcs, 50 pcs
098712	72x100 cm	10 pcs, 25 pcs, 50 pcs



*Available globally, for more information please contact info@sn-medical.com.

Alternative covers can be designed and manufactured per request.

C-Arm Covers

The C-Arm Cover is a functional, user-friendly, non-touch, fully sterile protection, and clip system, enabling full range of movement.

- Provides sterile cover of C-Arm.
- Clips allow effective, simple, and secure covering of the C-Arm system.
- Complete rotation and movement of the C-Arm, not compromising sterilization.
- Latex-free.
- Covers are made of high-grade clear polyethylene.

Siemens


For models: Siremobil / Multimobil / Arcadis / Cios Alpha.

3 Part Set: 1 clip system + 2 banded bags

Item No.	Description	Upper Cover	Lower Cover	Quantity
53ND5130	Sirmobil, Multimobil Series	75x90 cm	75x90 cm	8 pc/bx, 4 bx/cs
53ND5150	Arcadis	90x75 cm	75x75 cm	8 pc/bx, 4 bx/cs

2 Part Set: 1 clip system + 1 banded bag

Item No.	Description	Upper Cover	Lower Cover	Quantity
53ND5132	Sirmobil, Multimobil Series	75x90 cm	-	8 pc/bx, 4 bx/cs
53ND5152	Arcadis with Tapes	90x75 cm	-	8 pc/bx, 4 bx/cs



Philips

For models: BV 25, 26, 29, 200, 212, 300 / Libra / Endura / Pulsera / Veradius.

3 Part Set: 1 clip system + 2 banded bags

Item No.	Description	Upper Cover	Lower Cover	Quantity
53ND5200T	Pulsera with Tapes	90x75 cm	55x75 cm	10 pc/bx, 4 bx/cs
53ND5210	Libra with Tapes	75x90 cm	55x75 cm	10 pc/bx, 4 bx/cs
53ND5220	BV 25-300 / Endura / Libra	75x90 cm	55x75 cm	10 pc/bx, 4 bx/cs
53ND5220T	BV 25-300 / Endura / Libra with Tapes	65x130 cm	55x75 cm	10 pc/bx, 4 bx/cs
53ND5230	Veradius	65x130 cm	55x75 cm	10 pc/bx, 4 bx/cs



*Available globally, for more information please contact info@sn-medical.com.

Ziehm

For models: CB7D / 7000 / 8000 / Vision RFD, FD, 3D.

3 Part Set: 1 clip system + 2 banded bags

Item No.	Description	Upper Cover	Lower Cover	Quantity
53ND5320	CB7D	75x90 cm	75x75 cm	10 pc/bx, 4 bx/cs
53ND5340	7000, 8000 Series with Tapes	75x90 cm	75x90 cm	10 pc/bx, 4 bx/cs



2 Part Set: 1 clip system + 1 banded bag

Item No.	Description	Upper Cover	Lower Cover	Quantity
53ND5342	7000, 8000 Series with Tapes	75x90 cm	-	10 pc/bx, 4 bx/cs

GE/OEC

For models: D6 / D9 / 7000 / 8000 / 9000 / Mini 6000.

3 Part Set: 1 clip system + 2 banded bags

Item No.	Description	Upper Cover	Lower Cover	Quantity
53ND5550	7000/8000/9000 Series	90x75 cm	55x75 cm	8 pc/bx, 4 bx/cs
53ND5550T	7000/8000/9000 Series with Tapes	90x75 cm	55x75 cm	8 pc/bx, 4 bx/cs




*Available globally, for more information please contact info@sn-medical.com.

C-Arm Covers - Various

Suits all types of C-Arm equipment. Mini C-Arm Covers are suited for models: GE / Fluroscan / Orthoscan / Shimadzu. Functional, user-friendly, clip system permitting full range of movement.

Item No.	Description	Quantity
53ND5400	C-Arm Cover with Tapes 115x80 cm	20 pc/bx, 4 bx/cs
53ND5410	C-Arm Cover with Tapes 115x260 cm	15 pc/bx, 4 bx/cs
53ND5420	C-Arm Cover - Split - with Tapes 115x260 cm	15 pc/bx, 4 bx/cs
53ND5433	C-Arm Cover with Tapes 90x225 cm	20 pc/bx, 4 bx/cs
53ND22709	C-Arm Uni Protection Cover 90x225 cm	25 pc/bx, 3 bx/cs
53ND5152	C-Arm Cover - Arcadis - Set of 2 Parts	8 pc/bx, 4 bx/cs
53ND5153	C-Arm Cover - Universal	8 pc/bx, 4 bx/cs
53ND5420	C-Arm Cover - Split - with Tapes 115x260 cm	15 pc/bx, 4 bx/cs
53ND5433	C-Arm Long Cover with Tapes 90x250 cm	20 pc/bx, 4 bx/cs
53ND5132	C-Arm Cover - Sirmobile 2000 - Set of 2 Parts	8 pc/bx, 4 bx/cs
53ND5466	Mini C-Arm Cover with Rubber 138x160 cm	10 pc/bx, 4 bx/cs
53ND5460	Mini C-Arm Cover 138x170 cm	15 pc/bx, 4 bx/cs
53ND5555	C-Arm Cover - GE 9000 - Clips Only	8 pc/bx, 4 bx/cs



C-Bow Cover

A simple application, sterile cover attached with clips to the C-Bow of the C-Arm.

Item No.	Description	Upper Cover	Lower Cover	Quantity
53ND5600	Set of 3 Parts	75x75 cm	75x75 cm	20 pc/bx, 4 bx/cs
53ND5600T	Set of 3 Parts with Tapes	75x75 cm	55x75 cm	20 pc/bx, 4 bx/cs
53ND5602	Set of 2 Parts	65x130 cm	-	20 pc/bx, 4 bx/cs
53ND5603	Set of 3 Parts	80x100 cm	75x75 cm	20 pc/bx, 4 bx/cs

*Available globally, for more information please contact info@sn-medical.com.

Light Handle Covers

A single-use, sterile, light handle covers that fit easily and securely over operating room light handles.

Item No.	Description	Quantity
53NS01300	Light Handle Cover - Sterile	120 pc/bx, 4 bx/cs
51NS01300	Light Handle Cover - Non Sterile	Bulk
53NS01302	Light Handle Cover - Pack of 2 - Sterile	80 pc/bx, 4 bx/cs



Unicover®

A unique single-use, sterile, universal cover designed to provide a safe solution for all light handles and adapters.

Item No.	Description	Quantity
53NS0111001	Unicover - Sterile	60 pc/bx, 4 bx/cs
51NS0111001	Unicover - Non Sterile	Bulk



*Available globally, for more information please contact info@sn-medical.com.

Microscope Covers

High-Definition lens provides optimal non-distorted visualization.

- Velcro closures ensure a secure fit around oculars.
- Additional poly adhesive straps (green) for a tight fit around the microscope arm, securing excess cover material.
- Total protection and sterilization for operating room equipment.
- Latex-free.
- High-grade clear polyethylene.
- All microscope covers are standard-equipped with high-definition (HD) lenses.
- Glare Guard - the unique anti-glare angled lens eliminates unwanted reflection.





Zeiss & Moeller

Zeiss models: MD/ Pentero 800, 900 / Vario / Vario 700 / Sensera / Movena / Visu / Lumera i, 300, 700 / Rescan 700 / Pico.

Moeller models: Hi-R 700, 900, 1000 / 3-1000/ 20-1000 / Tocular / Allegra 90, 590, 900 / Universa 300 / Eibos 2.

Item No.	Description	Quantity
53ND76180	Ø6.5 - 52x180 cm - HD - Flat	25 pc/bx, 4 bx/cs
53ND86220	Ø6.5 - 115x220 cm - HD - Flat	15 pc/bx, 4 bx/cs
53ND86225	Ø6.5 - 115x220 cm - G.G. - Angle	7 pc/bx, 4 bx/cs
53ND86260	Ø6.5 - 115x260 cm - HD - Flat	12 pc/bx, 4 bx/cs
53ND86265	Ø6.5 - 115x260 cm - G.G. - Angle	7 pc/bx, 4 bx/cs
53ND86300	Ø6.5 - 115x300 cm - HD - Flat	12 pc/bx, 4 bx/cs
53ND96260	Ø6.5 - 137x260 cm - HD - Flat	10 pc/bx, 4 bx/cs
53ND96300	Ø6.5 - 137x300 cm - HD - Flat	10 pc/bx, 4 bx/cs
53ND96380	Ø6.5 - 137x380 cm - HD - Flat	8 pc/bx, 4 bx/cs
51ND96380	Ø6.5 - 137x380 cm - HD - Flat	Bulk
53ND96385	Ø6.5 - 137x380 cm - G.G. - Angle	5 pc/bx, 4 bx/cs
53ND87225	Ø7.0 - 115x220 cm - G.G. - Angle	7 pc/bx, 4 bx/cs



53ND86225

*Available globally, for more information please contact info@sn-medical.com.

Leica

Leica models: M220 F12 / M320 F12 / M400 E / M520, M525 OH4, F20, F40, F50 / M530 OH6 / M620 F20 OH3 / M651 OH1 / M720 OHS / M822 F40-F20/ M844 F40-F20 / MS-2, MS-3 / DI C500 / FL 400,800 / ULT 500.

Item No.	Description	Quantity
53ND87180	Ø7.0 - 115x180 cm - HD - Flat	15 pc/bx, 4 bx/cs
53ND87260	Ø7.0 - 115x260 cm - HD - Flat	12 pc/bx, 4 bx/cs
53ND87300	Ø7.0 - 115x300 cm - HD - Flat	12 pc/bx, 4 bx/cs
51ND87300	Ø7.0 - 115x300 cm - HD - Flat	Bulk
53ND97260	Ø7.0 - 137x260 cm - HD - Flat	10 pc/bx, 4 bx/cs
51ND97260	Ø7.0 - 137x260 cm - HD - Flat	Bulk
53ND97300	Ø7.0 - 137x300 cm - HD - Flat	10 pc/bx, 4 bx/cs
53ND97265	Ø7.0 - 137x260 cm - G.G. - Angle	6 pc/bx, 4 bx/cs
53ND97385	Ø7.0 - 137x380 cm - G.G. - Angle	5 pc/bx, 4 bx/cs
51ND97385	Ø7.0 - 137x380 cm - G.G. - Angle	Bulk



Dome Bags

The high-quality, sterile, single-use dome bag is designed to cover a wide range of equipment.

Item No.	Description	Quantity
53ND10050	50R 50x100 cm	40 pc/bx, 4 bx/cs
53ND10065	65R 65x130 cm	20 pc/bx, 4 bx/cs




*Available globally, for more information please contact info@sn-medical.com.

Banded Bags

The high-quality, sterile, single use Banded Bag provides a snug closure around the equipment.

Item No.	Description	Quantity
53ND17590	75x90 cm	50 pc/bx, 4 bx/cs
53ND17590T	75x90 cm with Tapes	50 pc/bx, 4 bx/cs
53ND18010	80x100 cm	50 pc/bx, 4 bx/cs
53ND18010T	80x100 cm with Tapes	50 pc/bx, 4 bx/cs
53ND18016	80x160 cm	40 pc/bx, 4 bx/cs
53ND11013	100x130 cm	30 pc/bx, 4 bx/cs



Camera Covers

The single-use Camera Cover is designed to provide a sterile interface between the unsterile operating room camera and the operating room staff.

Item No.	Description	Quantity
53ND2002	15x250 cm - Ring Type	30 pc/bx, 4 bx/cs
53ND02010	5x250 cm - Closed System	25 pc/bx, 4 bx/cs



Ultrasound Probe Covers

The extensive range Probe Covers guarantee maximum protection for the user, patient, and equipment.

Item No.	Description	Quantity
53ND2600	Standard - 2 Tapes - Telescopic Fold - 18x120 cm	320 pc/cs
53ND2705	2 Tapes - Telescopic Fold - 18x250 cm	60 pc/bx, 4 bx/cs
53A96090	20x80 cm	100 pc/cs

*Available globally, for more information please contact info@sn-medical.com.

Gowns

Basic Gown - Prevents infection from spreading and offers a sterile environment.


Standart Gown - Suitable for short medical procedures with minimal fluid involvement and for treating patients.

Reinforcement Gown - Offers a liquid-repellent fabric that is adapted for long procedures involving fluid-intensive environments, such as surgery and catheterization.

Fully Reinforcement Gown - A long-lasting gown, offering a high-grade liquid repellent, suitable for long surgeries involving a large amount of fluid.

Single Use Scrub Suits - Suitable for the use of the daily use of hospital staff, patients, and visitors. Protects from minimal fluid involvement.

Item No.	Description	Pkg.
53594520	Basic Gown - Sterile	1 pc/bag, 50 pc/cs
51594523	Basic Gown - Non Sterile	1 pc/bag, 50 pc/cs
53594400	Standard Gown (SMS) M	1 pc/bag, 25 pc/cs
53594410	Standard Gown (SMS) L	1 pc/bag, 25 pc/cs
53594420	Standard Gown (SMS) XL	1 pc/bag, 25 pc/cs
53594210	Reinf. Gown (SMS) L	1 pc/bag, 25 pc/cs
53594220	Reinf. Gown (SMS) XL	1 pc/bag, 25 pc/cs
53594230	Reinf. Gown (SMS) XXL	1 pc/bag, 25 pc/cs
53594325	Fully Reinf. Gown XL	25 pc/cs
51594502	Single Use Scrub Suit 2XL	50 pc/cs
51594503	Single Use Scrub Suit 3XL	50 pc/cs
51594505	Single Use Scrub Suit 4XL	50 pc/cs





*Available globally, for more information please contact info@sn-medical.com.

Masks

Sion Medical offers a variety of masks for several uses, such as masks for operating rooms, masks saturated with probiotics and other masks that filter out viruses and bacteria and provide partial protection from airborne diseases.

Item No.	Description	Pkg.
51578904	N95 Respirator Mask Type B	25 pc/bx, 15 bx/cs
51578903	N95 Respirator Mask Type A	25 pc/bx, 15 bx/cs
51578908-000	Level 3 Face Mask	50 pc/bx, 40 bx/cs
51578911	Probiotics Face Mask	50 pc/bx, 40 bx/cs
789898	Face Shield Mask	200 pc/cs
51592119	Protective Hood	50 pc/cs



*Available globally, for more information please contact info@sn-medical.com.