

Each set is available in different configurations and dimensions per the customer's request. Examples of sets appear on the following pages. Contact us for further specifications.



\*Available globally, for more information please contact [info@sn-medical.com](mailto:info@sn-medical.com).

## Ophthalmic

### Ophthalmic Pack

The ophthalmic set is suitable for eye surgeries such as cataract surgeries and pterygium surgeries. This set offers a variety of essential single-use items needed during most surgeries, such as pads and gauze, syringes and bowls, eye shields, and sort covers that are essential as a barrier between the sterile operation area and the contaminated area.

Description of component	Quantity
Gauze Ball G4	10
Sponge Forceps 19 cm	1
Tray 18x23 cm	1
Table Drape 75x100 cm	1
Applicator 6"	2
Cup 60ml	2
Polyethylene Bag 16x30 cm	1
Syringe 1LS	1
Needle 18GA	1
Syringe 10LL	1
Syringe 3LL	2
Syringe 5LL	1
Gauze Pad 8ply 7.5x7.5 cm	5
Section Irrigation Bag 30x50 cm	1
Eye Shield	1
Tray 13x17 cm	1
Eye Pad	1
Ophthalmic Drape 27.5x40 cm	1
Velcro Tube Holder 2.5x11.5 cm	1
Cup 60ml	2
Eye Applicator with Rubber	5
Cellulose Surgical Eye Spears	5
Image Cover 80x80 cm	2
Ophthalmic Drape 150x250 cm	1
Table Cover 100x140 cm	1
Standard Gown XL	2
Infusion Set Vented LS	1

\*Available globally, for more information please contact [info@sn-medical.com](mailto:info@sn-medical.com).

## Orthopedic

### Total Knee Replacement Pack

The set is intended for use in knee replacement surgeries. A complete set of covers is included in this set for all the equipment used in knee joint surgeries, such as a table cover and mayo cover, as well as patient covers to maintain a sterile environment during the operation, even when large amounts of fluid is involved.

Description of component	Quantity
Drape 180x200 cm with Adhesive	1
Mayo Cover 80x140 cm	2
Table Cover 137x220 cm	2
Velcro Tube Holder 2.5x11.5 cm	1
Foot Drape 170x200 cm	2
Stockinette S 30x75 cm	1
Split Sheet	1
Extremity Drape	1
O.P Tape 10x50 cm	1
Light Handle Cover	1
Reinf.Gown XL (SMS) + Towel	2

\*Available globally, for more information please contact [info@sn-medical.com](mailto:info@sn-medical.com).

## General Surgery

### Laparoscopy Pack

The Laparoscopy set is offered for minimally invasive pelvic cavity laparoscopy surgery. It offers a wide variety of medical equipment that is necessary for most procedures; including diathermia pencil, lap sponges, blades, camera cover, and other single-use covers such as gowns, table covers, mayo covers, and patient universal cover set that assures sterile operation area.

Description of component	Quantity
Lap Sponge 45x45 cm	10
Reinf. Gown XL + Towel	1
Reinf. Gown XXL + Towel	2
Suction Irrigation Bag 30x50 cm	1
Hand Towel 30x40 cm	2
Mayo Cover 60x190 cm	1
Table Cover 100x140 cm	1
Table Cover 152x250 cm	1
Camera Cover	1
Diathermy Pouch 2 Pocket 35x42 cm	1
Side Drape 75x100 cm	2
Foot Drape 170x200 cm	1
Head Drape 150x240 cm	1
Blade #11	1
Needle Counter	1
Diathermy Pencil 5 m	1
Suction Tube 350 cm	1
MIS HS-Spirit Medical Honey Dressing	1

\*Available globally, for more information please contact [info@sn-medical.com](mailto:info@sn-medical.com).

## Cardiology

### Cath Lab Pack

The Cardiovascular set is used for catheterization procedures. This set contains single-use products that are needed in every procedure, such as needles for the Selinger technique, Guidewire, syringes, needles, cups, blades, and a waste bowl used as a waste blood container. This set also includes covers that are essential as a barrier between the sterile operation area and the contaminated area.

Description of component	Quantity
Gauze Pad 12ply 10x10 cm	20
18GAX2.75 Introducer Needle	1
Guidewire 150 Length, 035 Diameter	1
Needle 18GA	1
Needle 21GA	1
Needle 25GA	1
Mini Scalpel #11	1
Blue Towel	4
Image Cover 80x80 cm	1
Protective Cover 80x100 cm	1
Radial Drape 280x350 cm	1
Table Cover 137x180 cm Without Reinf. Area	2
Syringe 3LL	2
Cup 60ml-Red	1
Cup 60ml-Green	1
Cup 60ml	2
Cup 100ml	2
Basin Guidewire 2500cc	1
Transparent Basin 250cc	1
Syringe 20LS	1
Syringe 20LL	1
Syringe 10LL	3
Waste Bowl	1
Sharps Holder	1
Blue Tray	1
Labels Set	1

\*Available globally, for more information please contact [info@sn-medical.com](mailto:info@sn-medical.com).

## Angio Pack

The Angio set is used for peripheral angiography catheterization procedures. This set offers single-use products needed for most Angio procedures, including covers for keeping a sterile procedure environment, and other single-use products as in the next list.

Description of component	Quantity
Gauze Pad 12 ply 10x10 cm	10
Gauze Ball G3	10
Sponge Forceps 19 cm	1
Infusion Set Vented LL 190 cm	1
Needle 25GA	1
Stopcock 4 Way	3
Blue Tower For Surgical Room 40x66 cm	4
Image Cover 80x80 cm	2
Key Board Cover 40x60 cm	1
Protective Cover 80x100 cm	1
Angio Drape 240x350 cm	1
Table Cover 137x180 cm	1
Needle 18GA	1
Needle 23G	1
Syringe 3LL	1
Scalpel #15	1
Needle Holder	1
Scalpel #11	1
Kocher Forceps	4
Extension 50 cm	1
Basin Guidewire 2500cc	2
Transparent Basin 250cc	2
Blue Basin 500cc	1
Syringe 20LL	3
Syringe 10LL	4
Sharps Holder	1
Kidney Basin 700cc	1

\*Available globally, for more information please contact [info@sn-medical.com](mailto:info@sn-medical.com).

## Neurosurgery Set

The Neurosurgery set is used in neurosurgery procedures. The set offers all covers needed for tables and other equipment needed for this operation. The patient drape has a transparent area for a better view of the surrounding area and a wide fluid pouch around the patient's head. This set contains sponges, laps, and other essential unique products such as suction tubes and skin staplers.

Description of component	Quantity
Gauze Pad 12 ply 10x10 cm XR	20
Lap Sponge 45x45 cm	10
Reinf. Gown XL + Towel	1
Skin Marker - Standard Tip	4
Tip Cleaner 4.5x4.5 cm	1
Light Handle Cover	6
Tray 19.5x22 cm	1
Blue Towel For Surgical Room 40x66 cm	8
Neuro Transparent Drape 220x240 cm	1
Mayo Cover 60x190 cm	1
Table Cover 100x140 cm	2
Table Cover 152x250 cm	1
Diathermy Pouch 2 Pocket 35x42 cm	2
Table Drape 75x100 cm	1
Line 14.5 cm	2
Needle 23g	2
Tommy Syringe 50cc	2
Syringe 10 LL	2
Syringe 20 LS	2
Cup 100cc	4
Scalpel #11	1
Blade #10	2
Skin Stapler	1
Needle Counter	1
Diathermy Pencil 5 m	1
Suction Tube 350 cm	3
Wide Sponge 6x6 cm	9
Venflon 20g	3

\*Available globally, for more information please contact [info@sn-medical.com](mailto:info@sn-medical.com).